

himac

APPLICATION

January 2003

Separation of cell organelle using easily prepared density gradient of Nycodenz®

CP-MX series preparative ultracentrifuge, P28S swinging bucket rotor

Experiment: Separation of rat liver cell organelle using the continuous density gradient without use of density gradient former

Cell organelle can be separated from a homogenate of liver etc. by the density gradient centrifugation using the ultracentrifuge. However, there was a tendency to avoid density gradient centrifugation because the conventional methods required the use of a density gradient former and the operation was complicated. This time, cell organelle were separated from a rat liver homogenate according to the procedure described in himac APPLICATION No. 108.

Continuous density gradient can be easily made according to the procedure below using Nycodenz® as the density gradient solution.

- (1) Dilute Nycodenz® to the specified concentration.
- (2) Freeze it at -20°C or -80°C .
- (3) Leave it at the ambient temperature to melt.

That is to say, you can obtain a continuous density gradient of Nycodenz® in the tube just by pouring the specified-concentration Nycodenz® in the tube, freezing it for preservation, taking it out of the freezer, and leaving it in the ambient temperature for melting. Separation of organelles in such a way is utilized as a part of proteome research and it seems to be important more and more in the future.

1. Equipment used

Centrifuge: CP-MX series preparative ultracentrifuge
 Rotor: P28S swinging bucket rotor
 Centrifugal tubes: 40PA tubes

2. Result of separation




Fig. 1 Result of electrophoresis (1-D SDS-PAGE) for each fraction after centrifugation (Amount of sample: $3 \mu\text{l}$)




Fig. 2. Western blotting result between fractionated samples and index proteins of each cell organelle

3. Separating conditions

Speed: 25,000 rpm (Average RCF: 82,200xg)

Time: 180 minutes

Temperature: 4°C

Sample: Rat liver homogenate

Amount of sample: 3.5 ml (Concentration: About 14mg/ml)

Density gradient solution: Nycodenz® (made by Nycomed Pharma AS (Oslo, Norway))

Initial concentration of Nycodenz®: 20%(w/v) (1.105g/ml)

Amount of Nycodenz®: 32 ml

Amount of fraction: 0.75 ml/Fraction

A fractionator is required for fraction after centrifugation. ALC-20 fractionator (made by Advantec Toyo Co., Ltd.) is reasonable and usable for fraction after centrifugation as well as Hitachi's DGF-U fractionator. Both fractionators are equipped with a level sensor and perform fractionation from the upper surface of the band in the centrifugal tube after centrifugation.

The above result was provided by Professor Kimie Murayama, Division of Biochemical Analysis, Central Laboratory of Medical Sciences, Juntendo University School of Medicine.

(References)

Kimie Murayama, Tsutomu Fujimura, Masataka Morita and Noriko Shindo, *Electrophoresis*, 2001, **22**, 2872-2880.

Nycodenz® is a trademark of Nycomed Pharma AS (Oslo, Norway).

HITACHI

Export offices

Hitachi High-Technologies Corporation

Head Office:

1-24-14, Nishishimbashi, Minato-ku, Tokyo, 105-8717 Japan
Tel: (81)3-3504-5461
Fax: (81)3-3504-7302

Manufacturer

Hitachi Koki Co., Ltd.

Scientific Instruments Division

1060, Takeda, Hitachinaka City, Ibaraki Pref., 312-8502 Japan
Tel: (81)29-276-7384
Fax: (81)29-276-7475

Nissei Sangyo (Singapore) Pte. Ltd.

Life Science Systems Dept.

10, Ang Mo Kio Street 65, #05-15 Techpoint, Singapore 569059
Tel: (65)481-2050
Fax: (65)481-8089

For the most current information, please access
<http://www.hitachi-koki.co.jp/himac/>